



## **HOME MEDICAL EQUIPMENT**

### **Education Information Checklist**

- **Where to place the Pap in the bedroom**
- **How to adjust and position the head gear and mask**
- **How to turn on the PAP unit**
- **How to activate the ramp function**
- **How to attach the mask, pillows, headgear and tubing**
- **How to change of clean the filters**
- **How to clean PAP unit and accessories**
- **How to disinfect the humidifier chamber**
- **Understand follow up program**
- **Refer to the manufacturers manual provided with your machine for further information**
- **Return Policy**

### **Safety**

- **Importance of grounding all electrical equipment**
- **How to obtain help if a medical emergency arises**
- **How to reach us for questions**

### **Cleaning Instructions**

- **For mask, tubing, and water chamber: Clean twice per week using Ivory/Dove dish soap and water. Allow to completely dry prior to use again**
- **Water chamber cleaning: Soak in 4-5 tablespoons of white vinegar and warm water for 30 minutes and then rinse**
- **Filter- Flick the dust off of this at least weekly**

### **Questions**

**Megan Mann, RRT 330-572-1011 x 4795**